

Hematoxylin and Eosin (H&E) Staining – Manual Protocol

(From Baylor College of Medicine)

Protocol for H&E staining:

- Place slides containing paraffin sections in a slide holder (glass or metal)
- Deparaffinize and rehydrate sections:
 - 3 x 3' Xylene (*blot excess xylene before going into ethanol*)
 - 3 x 3' 100% ethanol
 - 1 x 3' 95% ethanol
 - 1 x 3' 80% ethanol
 - 1 x 5' deionized H₂O
- While sections are in water, skim surface of hematoxylin with a Kimwipe to remove oxidized particles. Blot excess water from slide holder before going into hematoxylin.
- Hematoxylin staining:
 - 1 x 3' Hematoxylin
 - Rinse deionized water
 - 1 x 5' Tap water (*to allow stain to develop*)
 - Dip 8-12x (fast) Acid ethanol (*to destain*)
 - Rinse 2 x 1' Tap water
 - Rinse 1 x 2' Deionized water (*can leave overnight at this stage*)
 - (Blot excess water from slide holder before going into eosin)
- Eosin staining and dehydration:
 - 1 x 30 sec Eosin (*up to 45 seconds for an older batch of eosin*)
 - 3 x 5' 95% ethanol
 - 3 x 5' 100% ethanol (*blot excess ethanol before going into xylene*)
 - 3 x 15' Xylene
- Coverslip slides using Permount (xylene based).
- Place a drop of Permount on the slide using a glass rod, taking care to leave no bubbles.
- Angle the coverslip and let fall gently onto the slide. Allow the Permount to spread beneath the coverslip, covering all the tissue.
- Dry overnight in the hood.

Reagents for H&E staining:

- **Acid Ethanol:** 1 ml concentrated HCl + 400 ml 70% ethanol
- **Hematoxylin:** Poly Scientific (Bayshore, NY) #s212A
Harris hematoxylin with glacial acetic acid
- **Eosin:** Poly Scientific (Bayshore, NY) #s176
Eosin Phloxine stain, working
- **Permount:** Fisher Scientific #SP15-100
Histological mounting medium