St. Michael's

Inspired Care. Inspiring Science.

Immunofluorescence Staining Protocol

(from IHC world)

1. Preparation of Slides

A. Cell Lines

- Grow cultured cells on sterile glass cover slips or slides overnight at 37 º C
- Wash briefly with PBS
- Fix as desired. Possible procedures include:

10 minutes with 10% formalin in PBS (keep wet)

5 minutes with ice cold methanol, allow to air dry

5 minutes with ice cold acetone, allow to air dry

Wash in PBS

B. Frozen Sections

- Snap frozen fresh tissues in liquid nitrogen or isopentane pre-cooled in liquid nitrogen, embedded in OCT compound in cryomolds. Store frozen blocks at 80 °C.
- Cut 4-8 um thick cryostat sections and mount on superfrost plus slides or gelatin coated slides. Store slides at 80 °C until needed.
- Before staining, warm slides at room temperature for 30 minutes and fix in ice cold acetone for 5 minutes. Air dry for 30 minutes.
- Wash in PBS

C. Paraffin Sections

- Deparaffinize sections in xylene, 2x5min.
- Hydrate with 100% ethanol, 2x3min.
- Hydrate with 95% ethanol, 1min.
- Rinse in distilled water.
- Follow procedure for pretreatment as required.

2. Pretreatments of Tissue Sections

Antigenic determinants masked by formalin-fixation and paraffin-embedding often may be exposed by epitope umasking, enzymatic digestion or saponin, etc. Do not use this pretreatment with frozen sections or cultured cells that are not paraffin-embedded.

3. Procedure

- 1. Rinse sections in PBS-Tween 20 for 2x2min
- 2. **Serum Blocking**: incubate sections with normal serum block species same as secondary antibody, for 30 minutes to block non-specific binding of immunoglobulin. Note: since this protocol uses avidin-biotin detection system, avidin/biotin block may be needed based on tissue type. If you do, the avidin/biotin block should be done after normal serum block.
- 3. **Primary Antibody**: incubate sections with primary antibody at appropriate dilution in primary antibody dilution buffer for 1 hour at room temperature or overnight at 4 °C.
- 4. Rinse in PBS-Tween 20.
- 5. Secondary Antibody: incubate sections with biotinylated secondary antibody at appropriate dilution in PBS for 30 minutes at room temperature.
- 6. Rinse in PBS-Tween 20 for 3x2min.
- 7. **Detection:** incubate sections in FITC-Avidin D in PBS for 30 minutes at room temperature. Protecting slides from light starting from this step to the end by covering slides with aluminum foil or black box.
- 8. Rinse in PBS-Tween 20 for 3x2min.
- 9. Counterstain with PI or DAPI if desire.
- 10. Rinse in PBS-Tween 20.
- 11. Dehydrate through 95% ethanol for 2 min, 100% ethanol for 2x3min.
- 12. Coverslip with anti-fade mounting medium.